



## iPad®: LAMP Words for Life® Quick Reference Guide

For app support, visit [www.lampwflapp.com](http://www.lampwflapp.com)

For LAMP™ approach training, visit  
[www.aacandautism.com](http://www.aacandautism.com)

### Using Guided Access

1. On the iPad, select the **Settings** icon.
2. In the **General** settings, select **Accessibility**.
3. Scroll down and select **Guided Access**. Turn it **on**. You'll be prompted to set a passcode.
4. When in the LAMP WFL app, triple click the **Home** or **Side** button to activate Guided Access®. Triple click again to end the Guided Access session.

For additional information on using Guided Access, visit <https://support.apple.com/en-us/HT202612>

### Copying a Vocabulary File for Modification

A copy of an original file must be made before editing or using Vocabulary Builder. This keeps the original files in place in case they are needed at a later time.

1. Tap **Vocab** > **Choose New Vocab**.
2. Select the vocabulary file, then **Customize**.
3. Give the file a new name and description.
4. Tap **Save**.
5. The new file will open automatically.

### Setting up Multiple Languages

You may want to set up multiple languages to switch between while communicating.

- **Go to Vocabulary:** Use the Go to Vocabulary button action to link two custom vocabularies. This is a better option than Profiles, as you only need to program the button action on each vocab to link them back and forth. The voice automatically switches to the language of the vocabulary file.
- **Vocabulary Builder lists differences:** There will be a separate vocabulary builder list for each language, even if they are linked through a Go to Vocabulary button.

### Blocking & Restoring Access to the Editing Menus

*To password-protect the Editing menus:*

1. In the LAMP WFL app, choose **Menu** > **Settings**.
2. Under "Editing", turn **Use Menu Passcode** on.
3. Tap in the Passcode field, enter your password, and tap **Done**.

If you forget the password, enter **LAMPWFL** to overwrite the existing password.

*To remove editing from the Navigation Bar:*

1. Leave the TouchChat app, go to iOS **Settings**.
2. Select **Apps** on the left panel, then **LAMP WFL** on the right panel.
3. Disable **Allow Editing**.
4. Return to LAMP WFL. The **Vocab** and **Menu** text will disappear.

### Adding a Word to a Page

To add a word, an existing button may be edited or a new button created.

1. Navigate to the page where the new word will be added.
2. Select **Menu**.
3. Select **Edit Page**.
4. Select the button for the new word. If it is *blank*, select **Create New Button**; otherwise select **Edit This Button**.
5. Enter button preferences such as label, message, picture, button color, or actions.
6. When finished, select **Save**.
7. Select **Done** in the top left of the page to exit edit mode.

### Creating a New Page

1. From the app, select **Menu**.
2. Select **Edit Page**.
3. Select **Menu**.
4. Select **Add New Page**.
5. Name the page to represent its content.
6. If you are using a bilingual vocabulary, select the **Page Language**.
7. Select **Save**.
8. Add words to the blank page or select **Done** to add words at a later time.

## Creating a Button that Links to a New Page

---

Navigate to the page where you want to create the link.

1. Select **Menu**.
2. Select **Edit Page**.
3. Select the button that will link to the new page.
4. If the button is blank, select **Create New Button**; otherwise select **Edit This Button**.
5. Label the button (all capital letters). This automatically populates the **Button Message** field, so clear that text. Select an image and any Body Color (if needed).
6. Scroll down to the **Button Action** area and select **Add**.
7. Add a **Visit** action.
8. Select the page to link to or create a new page by selecting the “+” symbol.
9. Remove the **Speech Message** action by selecting **Edit** and choosing the “minus” button. Select **Delete**.
10. Select **Save**.
11. Select **Done**.

## Word Finder

---

Searching for a word is quick and easy with Word Finder. To use this feature:

1. From the app, select **Menu** or the **Keyboard**.
2. Select **Word Finder**.
3. Type the word you want to find.
4. Select **Find**. If the word is not stored, the **Find** button will remain grayed out and you will not be able to select it.
5. The window will identify which keys to press to activate the word.
6. Touch the icons to practice the motor plan.

## Changing the Voice

---

1. From the app, select **Menu**.
2. Select **Settings**.
3. In the **Speech** section, tap your voice to open the voice settings.
4. From the voice settings, tap **Voice** to open a list of available voices.
5. Select your preferred voice. Some voices require a free download by tapping the download icon.
6. Select **Back**, then **Done**.

## Using Vocabulary Builder

---

### Creating a List of Words to Teach

1. From the app, select **Menu**.
2. Select **Vocab Builder**.
3. Enable the **On/Off** button at the top left.
4. Tap in the Search field to open the keyboard.
5. Type a word (or partial word) you want to teach in the search window. Select the word you want from the list or tap “Search” on the keyboard.
6. You can display the Keyboard or Word Finder button by enabling Keyboard or Word Finder in the menu.
7. Select **Done** to return to the vocabulary.

## Managing Lists

At any time, a list may be saved, loaded, and shared.

1. From the app, select **Menu > Vocab Builder**

### Save a List

- a. Select **Save**
- b. Give the list a name in the **Exported Text** field.
- c. Select **Move**.

### Load a List

- a. Select **Load**.
- b. Select the preferred list to load from the Files menu.

### Share a List

- a. Select **Share** to share the currently loaded list.
- b. Select the method to share/back up the file.

2. Select **Done** to return to the vocabulary.

## Backing Up a Customized Vocabulary File

---

### Using AirDrop® or Email

1. On the app, select the **Vocab** button and select **Choose New Vocab**.
2. Select the **Menu** button and select **Import/Export Vocab**.
3. Select **Export and Share**.
4. Select the vocabulary you would like to export.
5. Select **Done** and select the method you would like to use to share.
  - a. **AirDrop**: Bluetooth® and Wi-Fi® must be active. You will see the nearby devices that have AirDrop on and can accept your file. Choose the desired device.
  - b. **Email**: Requires an email account set up on the iPad and access to Wi-Fi.

### Using iTunes or iShareWFL®

Please see the HELP section within the app.